

# A Single Calibration Method for Water AND Soil Samples Performing EPA Method 8260

Anne Jurek

## Introduction:

The United States Environmental Protection Agency (USEPA) Method 8260 is used to determine Volatile Organic Compounds (VOCs) in a variety of matrices. These matrices can vary from water to soil to sludge. In order to determine VOCs in the assorted matrices, the sampling and analysis system needs to be calibrated as closely to the matrix as possible. Therefore, USEPA Method 8260 requires a separate calibration for waters, soils and extractions. This application note will demonstrate a patent pending, automated water sampling mode using the soil sampling station of the autosampler thus eliminating the need to have separate calibrations for waters and soils.

## Discussion:

Environmental laboratories are always searching for ways to save time and run more samples per twelve hour increment. In order to save time, the bake time can be reduced, the purge time and flow can be changed, a dual sampling system can be run on one GC/MS, etc. However, one of the most basic ways to save time is by reducing the number of required curves and standards that need to be run. Unfortunately, the USEPA Method 8260 dictates the curves and sample experimental parameters need to be the same, and remain 5035 compliant, until now.


EST Analytical has designed a patent pending, automated water sampling mode that will allow environmental laboratories to run water samples in the soil mode. Thus, water and soil curves, standards, blanks and samples can use the same sampling parameters. (See steps 1, 2, and 3 below). This innovative sampling mode would eliminate the need to have separate curves, standards, etc. for soils and waters, thus saving laboratories time and money.



### Step 1:

Empty vials are placed in one tray while water standards or samples are placed in the second tray.

The system transports the empty vial to the soil sampling station and it is moved onto the sample needle.



### Step 2:

The arm moves over to the full water vial, the vial is pressurized with helium gas and the prescribed water volume is removed.

Internal standard (IS) is added by injecting IS into the sample as it is transported to the vial in the soil sampling station.



### Step 3:

The sample is heated and purged in the same manner as a soil sample, and the analytes are trapped on the analytical trap. Next, the analytes are desorbed onto the GC column.

## Experimental:

The sampling system used for this study was the EST Analytical Centurion WS autosampler and Encon Evolution concentrator. The concentrator was affixed with a Vocarb 3000 trap and connected to an Agilent 7890A GC and 5975 inert XL MS. The GC was configured with a Restek Rxi-624 Sil MS 30m x 0.25mm x 1.4µm column. The GC/MS experimental conditions are outlined in Table 1. The experiments were run using the extraction mode of the Centurion WS. Refer to Table 2 for the sampling method parameters.

GC/MS	Agilent 7890A/5975 inert XL
Inlet	Split/Splitless
Inlet Temp.	200°C
Inlet Head Pressure	12.153 psi
Mode	Split
Split Ratio	40:1
Column	Rxi-624Sil MS 30m x 0.25mm I.D. 1.4µm film thickness
Oven Temp. Program	45°C hold for 1 min., ramp 15°C/min to 220°C, hold for 1.3 min.
Column Flow Rate	1.0mL/min
Gas	Helium
Total Flow	44.00mL/min
Source Temp.	230°C
Quad Temp.	150°C
MS Transfer Line Temp.	180°C
Scan Range	m/z 35-265
Scans	3.12 scans/sec
Solvent Delay	0.7 min

**Table 1: GCMS Parameters**

Purge and Trap Concentrator	EST Encon Evolution
Trap Type	Vocarb 3000
Valve Oven Temp.	150°C
Transfer Line Temp.	150°C
Trap Temp.	35°C
Moisture Reduction Trap (MoRT) Temp.	39°C
Purge Time	11 min
Purge Flow	40mL/min
Dry Purge Temp.	ambient
Dry Purge Flow	40mL/min
Dry Purge Time	1.0 min
Desorb Pressure Control	On
Desorb Pressure	5psi
Desorb Time	0.5 min
Desorb Preheat Delay	5 sec.
Desorb Temp.	260°C
Moisture Reduction Trap (MoRT) Bake Temp.	230°C
Bake Temp	270°C
Sparge Vessel Bake Temp.	120°C
Bake Time	8
Bake Flow	40mL/min
Purge and Trap Auto-Sampler	EST Centurion WS
Sample Type	Water (Extraction Mode)
Sample Fill Mode	Syringe
Sample Volume	10mL
IS Volume	5µl
Needle Rinse Time	20 sec
Needle Sweep Time	20 sec
Concentrator Desorb Time	0.5 min
Syringe Rinse	On/12mL
Number of Syringe Rinses	2
Sparge Rinse Time	Off
Water Heater Temp.	85°C


**Table 2: Purge and Trap Water Parameters**

The USEPA Method 8260 standards were acquired from Restek. The linear range of the water extraction experiment was established by running a nine point calibration curve with a range of 0.5 to 200ppb. Method detection limits were also established for each compound by examining seven replicate standards of the low calibration point. Finally, precision and accuracy were determined by running seven replicate midpoint standards. Experimental Results are listed in Table 3.



Compound	Curve %RSD	Curve Ave. RF	MDL	Precision (%RSD 50ppb)	Accuracy Ave. 50ppb Recovery	Compound	Curve %RSD	Curve Ave. RF	MDL	Precision (%RSD 50ppb)	Accuracy Ave. 50ppb Recovery
Dichlorodifluoromethane	8.39	0.681	0.18	3.19	104.57	cis-1,3-Dichloropropene	10.18	0.824	0.08	1.05	108.02
Chloromethane	9.74	1.436	0.12	2.12	100.67	4-methyl-2-pentanone	3.11	0.621	0.08	1.98	104.54
Vinyl Chloride	8.59	1.529	0.09	1.41	112.97	Toluene-d8 SUR	9.15	2.026	0.09	1.10	108.73
Bromomethane	10.19	1.429	0.18	1.20	98.16	Toluene	8.6	1.383	0.09	0.94	110.91
Chloroethane	6.12	0.622	0.23	4.61	106.65	ethyl methacrylate	10.93	0.699	0.07	2.13	112.74
Trichlorofluoromethane	5.04	1.207	0.13	1.96	102.31	trans-1,3-Dichloropropene	11.07	0.803	0.12	0.90	111.66
diethyl ether	4.39	0.593	0.15	3.06	103.62	1,1,2-Trichloroethane	7.42	0.478	0.17	1.77	106.88
1,1,2-trichlorofluoroethane	3.9	0.731	0.19	1.68	102.41	Tetrachloroethene	11.34	0.641	0.22	0.96	125.02
1,1-Dichloroethene	4.26	0.680	0.10	2.02	102.43	1,3-Dichloropropane	6.91	0.833	0.10	1.45	106.40
Acetone	12.08	0.329	0.57	3.55	96.00	isopropyl acetate	10.84	0.141	0.23	1.85	109.14
Iodomethane	*0.999	0.540	0.10	3.31	98.65	butyl acetate	8.25	0.395	0.22	2.14	104.57
Carbon Disulfide	5.92	2.143	0.18	1.35	101.99	Dibromochloromethane	12.28	0.51	0.11	1.54	114.85
allyl chloride	6.59	0.997	0.18	1.87	103.38	2-Hexanone	4.51	0.469	0.25	2.21	103.87
Methylene Chloride	5.89	0.800	0.14	1.97	99.36	1,2-Dibromoethane	6.93	0.5	0.12	1.30	106.64
acetonitrile	11.85	0.092	0.37	5.85	93.77	Chlorobenzene	6.12	1.409	0.11	0.99	103.75
Tert Butyl Alcohol	9.98	0.144	1.09	3.46	102.11	1,1,1,2-Tetrachloroethane	12	0.492	0.15	1.62	113.41
MTBE	4.96	2.449	0.07	2.78	105.59	Ethylbenzene	8.85	2.499	0.12	1.10	110.51
cis-1,2-Dichloroethene	5.9	0.909	0.08	2.05	103.58	Xylene (m+p)	11.12	1.945	0.21	1.11	113.70
acrylonitrile	4.62	0.339	0.28	3.58	103.67	Styrene	13.39	1.559	0.09	1.27	115.55
Isopropylether	8.98	2.315	0.10	2.03	107.06	Xylene (o)	10.51	1.968	0.08	1.26	112.84
Vinyl acetate	9.98	2.712	0.04	2.67	95.55	n-amyl acetate	8.03	0.887	0.08	2.12	112.03
1,1-Dichloroethane	7.7	1.497	0.10	1.94	105.28	Bromoform	12.14	0.348	0.11	1.90	112.12
Ethyl Tert Butyl Ether (ETBE)	3.87	2.462	0.09	2.51	104.93	Isopropylbenzene	11.49	2.391	0.11	1.49	110.51
trans-1,2-Dichloroethene	4.26	0.680	0.10	2.02	102.43	cis-1,4-dichloro-2-butene	13.82	0.21	0.18	1.74	116.59
ethyl acetate	6.41	0.141	0.23	2.10	96.26	BFB SUR	8.61	1.101	0.15	1.95	99.76
2-Butanone	4.23	1.428	0.14	3.16	100.25	Bromobenzene	5.18	1.798	0.16	4.24	103.89
2,2-Dichloropropane	11.39	1.131	0.14	4.08	101.61	1,2,3-Trichloropropane	7.76	0.818	0.11	5.56	100.86
Bromochloromethane	5.99	0.503	0.09	1.80	102.76	1,1,2,2-Tetrachloroethane	6.76	0.99	0.07	2.55	100.40
propionitrile	9.02	0.156	0.21	2.28	104.45	n-Propylbenzene	7.36	4.695	0.09	1.83	105.83
methacrylonitrile	4.23	0.519	0.12	3.00	101.44	trans-1,4-dichloro-2-butene	6.24	0.355	0.23	2.30	103.67
THF	11.84	0.326	0.35	3.41	93.23	2-Chlorotoluene	4.95	0.941	0.12	1.76	100.84
Chloroform	4.79	1.455	0.18	1.57	104.32	4-Chlorotoluene	6.89	1.01	0.15	1.73	101.21
methyl acrylate	2.9	0.873	0.13	2.50	103.84	1,3,5-Trimethylbenzene	10.58	3.206	0.09	1.94	106.94
Dibromofluoromethane SUR	2.74	0.802	0.11	2.02	100.97	tert-Butylbenzene	11.97	2.827	0.11	1.98	103.40
1,1,1-Trichloroethane	4.04	1.256	0.08	1.36	105.53	sec-Butylbenzene	8.54	0.852	0.17	2.32	102.89
2-Chloroethylvinylether	8.86	0.342	0.18	2.89	120.49	1,2,4-Trimethylbenzene	9.52	3.286	0.13	1.88	106.27
Carbon Tetrachloride	12.96	1.038	0.09	1.84	109.80	nitrobenzene	*0.998	0.068	0.69	4.98	94.35
1,1-Dichloropropene	8.11	1.145	0.13	2.06	106.95	1,3-Dichlorobenzene	7.31	1.915	0.12	1.64	100.26
methyl acetate	6.01	0.859	0.10	2.50	102.39	1,4-Dichlorobenzene	10.59	1.999	0.18	2.04	96.53
Tert Amyl Methyl Ether (TAME)	5.67	2.370	0.07	2.80	107.09	Isopropyltoluene	13.34	3.486	0.11	1.91	107.66
Benzene	5.52	3.329	0.05	1.52	104.94	1,2-Dichlorobenzene	8.27	1.841	0.14	2.09	99.17
1,2-Dichloroethane	4.52	1.244	0.09	2.06	103.94	n-Butylbenzene	10.94	3.176	0.13	1.88	109.27
propyl acetate	6.23	0.884	0.13	1.63	105.03	1,2-Dibromo-3-chloropropane	11.97	0.219	0.35	3.12	107.17
Trichloroethene	6.86	0.328	0.10	0.57	103.48	1,2,4-Trichlorobenzene	10.62	1.275	0.13	2.12	95.00
1,2-Dichloropropane	10.27	0.544	0.11	1.25	105.25	Naphthalene	10.87	3.271	0.12	2.69	95.76
methyl methacrylate	9.17	0.474	0.08	1.82	111.39	Hexachlorobutadiene	7.03	0.529	0.13	2.77	100.38
Dibromomethane	6.64	0.351	0.11	1.40	103.25	1,2,3-Trichlorobenzene	11.84	1.125	0.10	2.35	97.25
Bromodichloromethane	9.86	0.745	0.08	1.07	109.68	Ave.	8.12	1.194	0.16	2.17	104.86
2-nitropropane	5.53	0.206	0.23	1.88	105.34						

Table 3: Experimental Results



**Figure 1: 50ppb Chromatogram of a Water Sample Run in Soil Mode**

**Conclusion:**

The patent pending sampling process of the Centurion WS proved to be a reliable and accurate sampling method for running water samples in the soil mode. The curve linearity and method detection limits both met USEPA Method 8260 requirements, and the precision and accuracy results were excellent. The water extraction option offered with the Centurion WS will save laboratories time and money because only one set of standards, curves, etc. is required for both water and soil samples.

**References:**

1. Volatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC/MS); United States Environmental Protection Agency Method 8260B, Revision 2, December 1996.

**Headquarters**

JSB International  
 Tramstraat 15  
 5611 CM Eindhoven  
 T +31 (0) 40 251 47 53  
 F +31 (0) 40 251 47 58

Zoex Europe  
 Tramstraat 15  
 5611 CM Eindhoven  
 T +31 (0) 40 257 39 72  
 F +31 (0) 40 251 47 58

**Sales and Service**

Netherlands  
 Apolloweg 2B  
 8239 DA Lelystad  
 T +31 (0) 320 87 00 18  
 F +31 (0) 320 87 00 19

Belgium  
 Grensstraat 7  
 Box 3 1831 Diegem  
 T +32 (0) 2 721 92 11  
 F +32 (0) 2 720 76 22

Germany  
 Max-Planck-Strasse 4  
 D-47475 Kamp-Lintfort  
 T +49 (0) 28 42 9280 799  
 F +49 (0) 28 42 9732 638

UK & Ireland  
 Cedar Court,  
 Grove Park Business Est.  
 White Waltham, Maidenhead  
 Berks, SL6 3LW  
 T +44 (0) 16 288 220 48  
 F +44 (0) 70 394 006 78

[info@go-jsb.com](mailto:info@go-jsb.com)  
[www.go-jsb.com](http://www.go-jsb.com)

With courtesy of

