


## Analysis of Gases Produces from the Pyrolysis of Biomass


Most biomass contains both cellulose and lignin, and when pyrolyzed using a capillary GC/MS produces a pyrogram like the one shown in Figure 1. This is a sample derived from pine wood, and the pyrolysis products include compounds from cellulose, like acetic acid and levoglucosan, and phenolics from the lignin such as methoxy phenol and vanillin (see Table 1 on the back).

Both cellulose and lignin produce fixed gases, including methane, carbon monoxide and water, which are difficult to analyze using capillary GC for several reasons. First, they may co-elute as the first peak, making specific identification difficult. Secondly, the masses for such small molecules may not be unique. Mass 28 could indicate nitrogen, carbon monoxide or ethylene. In addition, if a trapping system is being used, to facilitate slow heating, use of a reactor or heating in a reactive atmosphere, the fixed gases are likely to pass through the trap and be vented instead of being transferred to the GC.

In this analysis, a sample loop was placed after the sorbent trap, which collects the larger pyrolysis products, so that the fixed gases could be collected and transferred to a second GC with a packed column and thermal conductivity detector (TCD). Once the gases are inside the loop, a valve is rotated to flush them to the inlet of the packed column, producing the chromatogram shown in Figure 2. In this way, the effect on the production of these gases caused by changing parameters such as temperature or heating rate can be evaluated.



**Figure 1.** Py-GC/MS of pine wood biomass at 600°C.



**Figure 2.** Packed column-TCD analysis of gases from pyrolysis of pine biomass.

## ANALYTICAL PARAMETERS

### PYROLYSIS

Pyroprobe 5200

Interface: 325°C  
Ramp: 10°C/mS  
Temperature: 600°C  
Time: 15 seconds

### CAPILLARY GC/MS

Carrier: He  
Column: 5% Phenyl methyl silicone  
30 m x 0.25 mm  
Split: 50:1  
Program: 40°C for 2 minutes  
10°C/minute  
300°C for 5 minutes

### PACKED COLUMN GC/TCD

Column: Carboxen 1000, 60/80 mesh  
1/4" x 9'  
Carrier: He  
Flow: 30 ml/min  
Program: 50°C for 4 minutes  
30°C/minute  
230°C for 8 minutes

## TABLE 1

- 1 Acetic acid
- 2 Furfural
- 3 Phenol, 2-methoxy-
- 4 Levoglucosenone
- 5 Phenol, 2-methoxy-4-methyl-
- 6 2-Furancarboxaldehyde, 5-(hydroxymethyl)-
- 7 2-Methoxy-4-vinylphenol
- 8 Vanillin
- 9 Phenol, 2-methoxy-4-(1-propenyl)-, (E)-
- 10 Ethanone, 1-(4-hydroxy-3-methoxyphenyl)-
- 11 1,6-Anhydro- $\alpha$ -D-glucopyranose (levoglucosan)
- 12 Methyl-(2-hydroxy-3-ethoxy-benzyl)ether
- 13 4-Hydroxy-2-methoxycinnamaldehyde

FOR MORE INFORMATION  
CONCERNING THIS APPLICATION,  
WE RECOMMEND THE  
FOLLOWING READING:

A. A. Boateng, H. G. Jung and P. R. Adler, *Pyrolysis of energy crops including alfalfa stems reed canarygrass and eastern gamagrass*, Fuel 85 (2006) 2450.

© courtesy of CDS Analytical Inc.

### Headquarters

JSB International  
Tramstraat 15  
5611 CM Eindhoven  
T +31 (0) 40 251 47 53  
F +31 (0) 40 251 47 58

Zoex Europe  
Tramstraat 15  
5611 CM Eindhoven  
T +31 (0) 40 257 39 72  
F +31 (0) 40 251 47 58

### Sales and Service

Netherlands  
Apolloweg 2B  
8239 DA Lelystad  
T +31 (0) 320 87 00 18  
F +31 (0) 320 87 00 19

Belgium  
Grensstraat 7  
Box 3 1831 Diegem  
T +32 (0) 2 721 92 11  
F +32 (0) 2 720 76 22

Germany  
Max-Planck-Strasse 4  
D-47475 Kamp-Lintfort  
T +49 (0) 28 42 9280 799  
F +49 (0) 28 42 9732 638

UK & Ireland  
Cedar Court,  
Grove Park Business Est.  
White Waltham, Maidenhead  
Berks, SL6 3LW  
T +44 (0) 16 288 220 48  
F +44 (0) 70 394 006 78

[info@go-jsb.com](mailto:info@go-jsb.com)  
[www.go-jsb.com](http://www.go-jsb.com)

